

Product datasheet for **AM08017HR-N**

Mouse IgG2a (gamma 2a chain specific) Rat Monoclonal Antibody [Clone ID: SB84a]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a (gamma 2a chain specific) Rat Monoclonal Antibody [Clone ID: SB84a]
Clone Name:	SB84a
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000-1/2000. Also used as a second step reagent for mouse IgG2a monoclonal antibodies.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)