

Product datasheet for **AM08017RP-S**

Mouse IgG2a (gamma 2a chain specific) Rat Monoclonal Antibody [Clone ID: SB84a]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG2a (gamma 2a chain specific) Rat Monoclonal Antibody [Clone ID: SB84a]
Clone Name:	SB84a
Applications:	ELISA, FC, IF, WB
Recommended Dilution:	Flow Cytometry: $\leq 0.2 \mu\text{g}/10^6$ cells. Identification and enumeration of IgG2a+. Immunofluorescence Microscopy: Identification and enumeration of IgG2a+. ELISA/FLISA. Also used as a second step reagent for mouse IgG2a monoclonal antibodies.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)