

## Product datasheet for **AM08018BT-N**

### Mouse IgA (alpha chain specific) Rat Monoclonal Antibody [Clone ID: 11-44-2]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA (alpha chain specific) Rat Monoclonal Antibody [Clone ID: 11-44-2]
Clone Name:	11-44-2
Applications:	ELISA, FC, IF
Recommended Dilution:	<b>Flow Cytometry:</b> < / = 1 µg/10e6 cells. Identification and enumeration of IgA+. <b>Immunofluorescence Microscopy:</b> Identification and enumeration of IgA+. <b>ELISA.</b>
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)