

Product datasheet for **AM08018FC-N**


Mouse IgA (alpha chain specific) Rat Monoclonal Antibody [Clone ID: 11-44-2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA (alpha chain specific) Rat Monoclonal Antibody [Clone ID: 11-44-2]
Clone Name:	11-44-2
Applications:	FC, IF
Recommended Dilution:	FLISA: 1/200-1/400. Flow Cytometry: $\leq 1 \mu\text{g}/10^6$ cells. Identification and enumeration of IgA+ cells. Immunofluorescence Microscopy: Identification and enumeration of IgA+ cells.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)

Product images:

FLISA plate was coated with purified mouse IgM, IgG, and IgA. Immunoglobulins were detected with serially diluted Rat Anti-Mouse IgA-FITC (Cat. AM08018FC-N).