

Product datasheet for **AM08018PU-N**

Mouse IgA (alpha chain specific) Rat Monoclonal Antibody [Clone ID: 11-44-2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA (alpha chain specific) Rat Monoclonal Antibody [Clone ID: 11-44-2]
Clone Name:	11-44-2
Applications:	ELISA, FC, IF
Recommended Dilution:	ELISA. Flow Cytometry: Identification and enumeration of IgA+. Immunofluorescence Microscopy: Identification and enumeration of IgA+.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)