

Product datasheet for **AM08018RP-S**

Mouse IgA (alpha chain specific) Rat Monoclonal Antibody [Clone ID: 11-44-2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgA (alpha chain specific) Rat Monoclonal Antibody [Clone ID: 11-44-2]
Clone Name:	11-44-2
Applications:	FC, IF
Recommended Dilution:	Flow Cytometry: < / = 0.2 µg/10e6 cells. Identification and enumeration of IgA+ cells. Immunofluorescence Microscopy: Identification and enumeration of IgA+ cells.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)