

Product datasheet for **AM08019AP-N**

Mouse Kappa Light Chain Rat Monoclonal Antibody [Clone ID: 187.1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	187.1
Applications:	ELISA
Recommended Dilution:	ELISA: 1/500-1/1000.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody reacts with Mouse kappa light chains. Does not react with other immunoglobulin isotypes or lambda light chains.
Formulation:	50mM Tris/1mM MgCl ₂ /50% Glycerol, pH 8.0 containing 0.09% Sodium Azide as preservative. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)