

## Product datasheet for **AM08019HR-N**

### Mouse Kappa Light Chain Rat Monoclonal Antibody [Clone ID: 187.1]

#### Product data:

|                       |   |
|-----------------------|---|
| Product Type:         | Primary Antibodies  |
| Clone Name:           | 187.1   |
| Applications:         | ELISA   |
| Recommended Dilution: | <b>ELISA.</b> 1/4000-1/8000.  |
| Reactivity:           | Mouse   |
| Host:                 | Rat   |
| Isotype:              | IgG1  |
| Clonality:            | Monoclonal  |
| Specificity:          | This antibody reacts with Mouse kappa light chains.<br>Does not react with other immunoglobulin isotypes or lambda light chains.  |
| Formulation:          | 50% Glycerol/50% PBS, pH 7.4 without preservatives.<br>Label: HRP<br>State: Liquid purified Ig fraction.<br>Label: Horseradish Peroxidase   |
| Conjugation:          | HRP   |
| Storage:              | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>This product is photosensitive and should be protected from light.<br>Avoid repeated freezing and thawing. |
| Stability:            | Shelf life: one year from despatch.   |



[View online »](#)