

## Product datasheet for **AM08019PU-N**

### Mouse Kappa Light Chain Rat Monoclonal Antibody [Clone ID: 187.1]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	187.1
Applications:	ELISA, FC, IF, IP
Recommended Dilution:	<b>Flow Cytometry:</b> Quality tested (Ref.4-15) <b>ELISA:</b> Quality tested (Ref.2,3). <b>FLISA:</b> Quality tested. <b>Immunohistochemistry on Frozen Sections:</b> (Ref.16). <b>Immunoprecipitation</b> (Ref.19). <b>Purification:</b> (Ref.20,21). <b>Costim:</b> Reported.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody reacts with Mouse kappa light chains. Does not react with other immunoglobulin isotypes or lambda light chains.
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.



[View online »](#)