

Product datasheet for **AM08019RP-S**

Mouse Kappa Light Chain Rat Monoclonal Antibody [Clone ID: 187.1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	187.1
Applications:	FC, IF
Recommended Dilution:	Identification and enumeration of kappa cells by Flow Cytometry: < = 0.2 µg/10e6 cells. Immunofluorescence.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	MPC-11 immunoglobulin
Specificity:	This antibody reacts with Mouse kappa light chains. Does not react with other immunoglobulin isotypes or lambda light chains.
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)