

Product datasheet for **AM08020FC-N**

Mouse Kappa Light Chain Rat Monoclonal Antibody [Clone ID: H139-52.1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	H139-52.1
Applications:	ELISA, FC, IF
Recommended Dilution:	Flow Cytometry: < / = 1 µg/10e6 cells. Identification and enumeration of kappa+ cells. Immunofluorescence Microscopy: Identification and enumeration of kappa+ cells.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody reacts with Mouse kappa light chains.
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)