

## Product datasheet for **AM08021HR-N**

### Mouse lambda light chain Rat Monoclonal Antibody [Clone ID: JC5-1]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	JC5-1
Applications:	ELISA
Recommended Dilution:	<b>ELISA.</b> 1/1000-1/2000.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2b
Clonality:	Monoclonal
Specificity:	This antibody reacts with mouse Lambda 1 and Lambda 2 light chains. It does not react with other immunoglobulin isotypes or Kappa light chains.
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)