

Product datasheet for **AM08021PU-N**

Mouse lambda light chain Rat Monoclonal Antibody [Clone ID: JC5-1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	JC5-1
Applications:	ELISA, FC, IF
Recommended Dilution:	Flow Cytometry: Identification and enumeration of Lambda+ cells. Immunofluorescence Microscopy: Identification and enumeration of Lambda+ cells. ELISA.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2b
Clonality:	Monoclonal
Specificity:	This antibody reacts with mouse Lambda 1 and Lambda 2 light chains. It does not react with other immunoglobulin isotypes or Kappa light chains.
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)