

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Product datasheet for AM08023RC7-N

## Cd8a Rat Monoclonal Antibody [Clone ID: 53-6.7]

## **Product data:**

Product Type:	Primary Antibodies
Clone Name:	53-6.7
Applications:	FC, IHC, IP
Recommended Dilution:	Flow Cytometry: < / = 0.2 μg/10e6 cells. (Ref.8,9) Immunohistochemistry (Acetone-Fixed Frozen Sections) (See Ref. 10) Immunoprecipitation. (Ref.8) Fractionation of CD8a+ T cells. (Ref.10)
Reactivity:	Mouse
Host:	Rat
lsotype:	lgG2a
Clonality:	Monoclonal
Specificity:	This antibody recognises the alpha chain of the Murine CD8 hetrerodimer (Mr. 32-34 kDa).
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE-Cy7 State: Liquid purified Ig fraction. Label: R-Phycoerythrin-Cyanine 7
Concentration:	lot specific
Conjugation:	PE-Cy7
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD8 antigen, alpha chain
Database Link:	Entrez Gene 12525 Mouse P01731



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Cd8a Rat Monoclonal Antibody [Clone ID: 53-6.7] – AM08023RC7-N
Background:	In the mouse, CD8 exists in two forms: (i) a CD8 heterodimer composed of an alpha chain (CD8a/Lyt-2) and a beta chain (CD8b/Lyt-3); and (ii) a homodimer of two alpha chains. The heterodimer is found on the surface of essentially all thymocytes and the "suppressor/cytotoxic" subpopulation of mature T lymphocytes. Subsets of intestinal intraepithelial lymphocytes express CD8a without CD8b. It has been suggested that CD8+b-T cells mature extrathymically, while development of the CD8a+b+ population of T cells is thymus-dependent. CD8 acts as a coreceptor with MHC Class I-restricted T cell receptors in antigen recognition and positive selection of MHC class I-restricted CD8+ T cells. In vivo and in vitro treatment with the 53-6.7 monoclonal antibody effectively depletes CD8a+ cells. The 53-6.7 monoclonal antibody also blocks allogeneic help specific for class I MHC antigens and T cell responses to IL-2. (Ref. 1-8)
Synonyms:	CD8 alpha chain, CD8A, MAL

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US