

Product datasheet for **AM08025LE-N**

Cd2 Rat Monoclonal Antibody [Clone ID: 12-15]

Product data:

Product Type:	Primary Antibodies
Clone Name:	12-15
Applications:	FC, FN, IP
Recommended Dilution:	Flow Cytometry < / = 1 µg/10e6 cells. (Ref.1-3) Immunoprecipitation. (Ref.1-2) Functional Assay: <i>In vitro</i> Inhibition of CTL-mediated lysis. (Ref.1) Induction of T cell unresponsiveness <i>in vivo</i> . (Ref.4)
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody is specific to CD2/LFA-2 (Mr 45-58 kDa)
Formulation:	PBS without preservatives. State: Low Endotoxin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD2 antigen
Database Link:	Entrez Gene 12481 Mouse P08920



[View online »](#)

Background:

CD2, also known as leukocyte function-associated antigen-2 (LFA-2), is a member of the immunoglobulin superfamily of cell adhesion molecules. (Ref.1) It is expressed on the surface of >95% of thymocytes, and essentially all B cells, T cells and NK cells. (Ref.1-3) In the mouse, CD2 is the receptor for CD48. (Ref.4) CD2 is thought to play a major role in mediating T-cell adhesion, T-cell recognition/activation, T-cell signaling and intrathymic differentiation. (Ref.5-10)

Synonyms:

SRBC, Erythrocyte receptor, LFA-2, LFA-3 receptor, Rosette receptor