

Product datasheet for AM08027PU-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Collagen IV (COL4A1) Mouse Monoclonal Antibody [Clone ID: 2F11]

Product data:

Product Type: Primary Antibodies

Clone Name: 2F11

Applications: ELISA, IHC

Recommended Dilution: ELISA: 1/2,000-1/8,000.

Immunohistochemistry: < / = 2 µg/ml.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Specificity: This antibody reacts with Human Type IV Collagen.

Formulation: 100 mM Borate Buffered Saline, pH 8.2.

No preservatives or amine-containing buffer salts added.

State: Purified

State: Liquid purified Ig fraction.

Concentration: lot specific

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: collagen type IV alpha 1 chain

Database Link: Entrez Gene 1282 Human

P02462

Background: Collagen IV is a major constituent of the basement membranes along with laminins,

proteoglycans and enactins. It is a multimeric protein composed of 3 alpha subunits. These subunits are encoded by 6 different genes, alpha 1 through alpha 6, each of which can form a triple helix structure with 2 other subunits to form type IV collagen. It can form insoluble fibers with high tensile strength. Collagen IV is useful in detecting the loss of parts of

basement membranes in carcinomas.





Collagen IV (COL4A1) Mouse Monoclonal Antibody [Clone ID: 2F11] - AM08027PU-N

Synonyms: COL4A1

Protein Pathways: ECM-receptor interaction, Focal adhesion, Pathways in cancer, Small cell lung cancer