

Product datasheet for **AM08027RP-N**

Collagen IV (COL4A1) Mouse Monoclonal Antibody [Clone ID: 2F11]

Product data:

Product Type:	Primary Antibodies
Clone Name:	2F11
Applications:	ELISA, IHC
Recommended Dilution:	ELISA: 1/2,000-1/8,000. Immunohistochemistry: < / = 2 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody reacts with Human Type IV Collagen.
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	collagen type IV alpha 1 chain
Database Link:	Entrez Gene 1282 Human P02462



[View online »](#)

Background:	Collagen IV is a major constituent of the basement membranes along with laminins, proteoglycans and enactins. It is a multimeric protein composed of 3 alpha subunits. These subunits are encoded by 6 different genes, alpha 1 through alpha 6, each of which can form a triple helix structure with 2 other subunits to form type IV collagen. It can form insoluble fibers with high tensile strength. Collagen IV is useful in detecting the loss of parts of basement membranes in carcinomas.
Synonyms:	COL4A1
Protein Pathways:	ECM-receptor interaction, Focal adhesion, Pathways in cancer, Small cell lung cancer