

## Product datasheet for AM08032AC7-N

## Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

**OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

## Cd19 Mouse Monoclonal Antibody [Clone ID: MB19-1]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: MB19-1

Applications: FC

**Recommended Dilution:** Flow Cytometry: < / = 0.2 µg/10e6 cells. (Ref.1-3)

Reactivity: Mouse
Host: Mouse
Isotype: IgA

Clonality: Monoclonal

Immunogen: CD19+ mouse pre-B cell line, 300.19. (Ref.1)

**Specificity:** This antibody is specific to CD19 B cell differentiation antigen (Mr. 95 kDa). **Formulation:** PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.

Label: APC-Cy7

State: Liquid purified Ig fraction. Label: Allophycocyanin-Cyanine 7

**Concentration:** lot specific **Conjugation:** APC-Cy7

Storage: Store the antibody undiluted at 2-8°C.

**DO NOT FREEZE!** 

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

Gene Name: CD19 antigen

**Database Link:** Entrez Gene 12478 Mouse

P25918





## Cd19 Mouse Monoclonal Antibody [Clone ID: MB19-1] - AM08032AC7-N

Background:

CD19 is a monomeric transmembrane glycoprotein expressed at relatively constant levels throughout B cell development from early pro-B/pre-B cells (i.e., B220+/CD43+/HSA+) through fully differentiated B cell stages. (Ref.1-7) Terminally differentiated plasma cells do not express CD19.1 In humans, the CD19 molecule on the surface of mature B cells associates with CD21 (CR-2) and CD81 (TAPA-1), and this multimolecular complex synergizes with surface immunoglobulin to provide signal transduction and promote cellular activation. (Ref.3,4) All splenic and peritoneal IgM+ cells of both B-1 and B-2 lineages are CD19+, with B-1 cells expressing higher levels of CD19 than B-2 cells in these sites. (Ref.1,7) Recent studies with CD19-deficient mice have suggested that this molecule may not be required for normal generation and maturation of B cells in the bone marrow. (Ref.7)

Synonyms: Leu-12, B-cell marker