

## Product datasheet for **AM08037AC-N**

### Cd28 Hamster Monoclonal Antibody [Clone ID: 37.51]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	37.51
Applications:	FC
Recommended Dilution:	<b>Flow Cytometry:</b> < / = 0.2 µg/10e6 cells. (Ref.2-5)
Reactivity:	Mouse
Host:	Hamster
Isotype:	IgG
Clonality:	Monoclonal
Specificity:	This antibody is specific to CD28 costimulatory receptor (Mr. 40 kDa)
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: APC State: Liquid purified Ig fraction. Label: Allophycocyanin
Concentration:	lot specific
Conjugation:	APC
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD28 antigen
Database Link:	<a href="#">Entrez Gene 12487 Mouse P31041</a>



[View online »](#)

- Background:** CD28 is a type I disulfide-linked homodimer that is constitutively expressed on most thymocytes, at low density on nearly all CD4+ and CD8+ peripheral T lymphocytes, and at very low levels on NK cells. Its expression is upregulated upon T-cell activation. (Ref.1-4) CD28 is a ligand for CD80/B7-1 and CD86/B7-2 on B cells and other antigen presenting cells, and plays an important role in the interaction between T cells and B cells. CD28 is a co-stimulatory receptor involved in many, but not all, T-cell independent immune responses. (Ref.5-8)
- Synonyms:** TP44