

Product datasheet for AM08041LE-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Cd38 Rat Monoclonal Antibody [Clone ID: NIMR-5]

Product data:

Product Type: Primary Antibodies

Clone Name: NIMR-5

Applications: FC, FN, IHC, IP

Recommended Dilution: Flow Cytometry. (Ref. 1-8)

Immunoprecipitation. (Ref.1)

Immunohistochemistry on Frozen Sections. (Ref.7)

Functional Assay: Induction of B lymphocyte proliferation. (Ref.3-8)

Reactivity: Mouse
Host: Rat
Isotype: IgG2a

Clonality: Monoclonal

Specificity: This antibody is specific to CD38 (Mr. 42-46 kDa).

Monoclonal antibodies to CD38 have been shown to induce B- and T-cell proliferation, protect

B cells from apoptosis, and inhibit B lymphopoiesis. (Ref.1-8)

Formulation: PBS containing no preservatives.

State: Low Endotoxin

State: Liquid purified Ig fraction.

Concentration: lot specific

Conjugation: Unconjugated

Storage: Store the antibody undiluted (in aliquots) at -20°C.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: CD38 antigen

Database Link: Entrez Gene 12494 Mouse

P56528





Cd38 Rat Monoclonal Antibody [Clone ID: NIMR-5] - AM08041LE-N

Background:

Murine CD38, a type II transmembrane glycoprotein, is a bifunctional ectoenzyme capable of catabolizing Nicotinamide adenine dinucleotide (NAD+) to cyclic ADP-ribose (cADPR), and then hydrolyzing cADPR to adenosine diphosphoribose (ADPR). It is expressed at high levels on the surface of peripheral B-lineage cells, and at low density on germinal center B cells from unimmunized mice. It has also been reported to be expressed at moderate levels on NK cells, a proportion of peripheral T cells, and a subpopulation of thymocytes which are mostly TCRalphabeta+, CD4-, CD8-. Murine CD38 is also expressed by all Mac-1+ macrophages in the peritoneal cavities of unimmunized mice, but not by unstimulated bone-marrow-derived macrophages.

Synonyms:

cADPr hydrolase 1, T10