

Product datasheet for AM08044RC5-N

CD45 / LCA Rat Monoclonal Antibody [Clone ID: 13/2.3]

Product data:

Product Type: Primary Antibodies Clone Name: 13/2.3IHC **Applications: Recommended Dilution: Flow Cytometry:** < / = 0.2 µg/10e6 cells. (Ref.4,5) **Reactivity:** Mouse Host: Rat Isotype: lgG2b **Clonality:** Monoclonal This antibody recognizes a framework epitope present on all CD45 isoforms. (Ref.3-6) Specificity: PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Formulation: Label: PE-Cy5.5 State: Liquid purified lg fraction. Label: R-Phycoerythrin-Cyanine 5.5 Concentration: lot specific **Conjugation:** PE-Cy5.5 Store the antibody undiluted at 2-8°C. Storage: **DO NOT FREEZE!** This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Stability: Shelf life: one year from despatch. Gene Name: protein tyrosine phosphatase, receptor type, C Entrez Gene 19264 Mouse Database Link: P06800



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

	CD45 / LCA Rat Monoclonal Antibody [Clone ID: I3/2.3] – AM08044RC5-N
Background:	CD45 is a family of single chain transmembraneous glycoproteins consisting of at least four isoforms (220, 205, 190, 180 kDa) which share a common large intracellular domain. Their extracellular domains are heavily glycosylated. The different isoforms are produced by alternative messenger RNA splicing of three exons of a single gene on chromosome 1. CD45 is expressed on cells of the human hematopoietic lineage (including hematopoietic stem cells) with the exception of mature red cells. It is not detected on differentiated cells of other tissues. It is likely that CD45 plays an important role in signal transduction, inhibition or upregulation of various immunological functions. Antibodies recognising a common epitope on all of the isoforms are termed CD45 while those recognising only individual isoforms are termed CD45RA or CD45RO etc.
Synonyms:	PTPRC, Leukocyte common antigen, L-CA, T200

%// _

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US