

Product datasheet for AM08054PU-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Thy1 Rat Monoclonal Antibody [Clone ID: G7]

Product data:

Product Type: Primary Antibodies

Clone Name: G7

Applications: FC, FN, IHC

Recommended Dilution: Flow Cytometry. (Ref.1,10)

Immunohistochemistry (Acetone-Fixed, Frozen Sections). (Ref. 10)

In vitro **T Cell Activation.** (Ref.1,5,6,9)

Induction of apoptosis.

Reactivity: Mouse
Host: Rat
Isotype: IgG2c

Clonality: Monoclonal

Specificity: This antibody recognises Thy-1.1 and Thy-1.2 alloantigens (Thy-1 epitope region A).

MAb G7 stimulates T-cell proliferation and IL-2 secretion, via signalling through the T-cell

receptor/CD3 complex. (Ref.1,5-7).

It has also been reported to promote apoptosis of thymocytes and CTL clones (Ref.4,7) and to

mediate adhesion of thymocytes to thymic stroma. (Ref.8)

Formulation: 100 mM Borate Buffered Saline, pH 8.2.

No preservatives or amine-containing buffer salts added.

State: Purified

State: Liquid purified Ig fraction.

Concentration: lot specific
Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: thymus cell antigen 1, theta

Database Link: Entrez Gene 21838 Mouse

P01831





Thy1 Rat Monoclonal Antibody [Clone ID: G7] - AM08054PU-N

Background: CD90/Thy-1, a GPI-anchored molecule and one of the smallest members of the

immunoglobulin superfamily of cell surface receptors, consists of a single V-set domain. (Ref.1-4) It is expressed on thymocytes, peripheral T lymphocytes, some intraepithelial T

lymphocytes and neurons of all mouse strains. (Ref.1,2)

Synonyms: Thy-1, THY1, CDw90