

Product datasheet for **AM08076RP-N**

MHC Class II I-Ab Mouse Monoclonal Antibody [Clone ID: 25-9-3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	25-9-3
Applications:	FC
Recommended Dilution:	Flow Cytometry: < / = 0.2 µg/10e6 cells.
Reactivity:	Mouse
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Immunogen:	C3H.SW splenocytes. (Ref.1)
Specificity:	This antibody reacts specifically with the I-Ab haplotype of MHC Class II molecules. Class II antigens are predominantly expressed on antigen-presenting cells including B lymphocytes, macrophages, dendritic cells and certain epithelial cells.
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P18468



[View online »](#)

Background: MHC Class II antigens are involved in antigen presentation by Antigen Presenting Cells (APCs) to CD4+ T cells. They are constitutively expressed on APCs including B cells, macrophages, monocytes and dendritic cells, and are inducible on a number of other cells (endothelium and epithelial cells) by interferon gamma. The Mouse H2 Ab locus is orthologous to human DQB, which varies from typical class II genes in that both the alpha and beta chains are polymorphic.

Synonyms: H2-Eb1, H-2 class II histocompatibility antigen I-A beta chain