

Product datasheet for **AM08077RP-N**

MHC Class II I-Ek Mouse Monoclonal Antibody [Clone ID: 17-3-3S]

Product data:

Product Type:	Primary Antibodies
Clone Name:	17-3-3S
Applications:	FC
Recommended Dilution:	Flow Cytometry: < / = 0.2 µg/10e6 cells. (Ref.1)
Reactivity:	Mouse
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	C3H mouse skin graft and splenocytes. (Ref.1)
Specificity:	This antibody is An epitope on mouse I-Ek MHC class II alloantigen. It reacts with the I-Ek class II alloantigen on cells from mice of the H-2r haplotype. Reactivity with H-2d and H-2p haplotypes has not been observed. (Ref.1) Cells from mice of the H-2b, H-2f, H-2g, H-2q, and H-2s haplotypes do not express I-E antigen. (Ref.2)
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)