

Product datasheet for **AM08081PU-S**

MHC Class I H2 Kd/Dd Mouse Monoclonal Antibody [Clone ID: 34-1-2S]

Product data:

Product Type:	Primary Antibodies
Clone Name:	34-1-2S
Applications:	CT, FC, IHC
Recommended Dilution:	Immunohistochemistry on frozen sections. Cytotoxicity Assay. Flow Cytometry.
Reactivity:	Mouse
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	BDF splenocytes <u>Donor:</u> C3H splenocytes <u>Fusion Partner:</u> Sp2/0-Ag14
Specificity:	This anti H-2KdDd mAb reacts with both H-2Kd and H-2Dd products. The antibody also cross reacts with Kbsrpq. This K-D cross reaction indicates the presence of shared specificities between the two separate H-2 regions.
Formulation:	State: Purified State: Liquid purified Ig fraction
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)