

## Product datasheet for AM08081RP-N

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## MHC Class I H2 Kd/Dd Mouse Monoclonal Antibody [Clone ID: 34-1-2S]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: 34-1-2S Applications: FC, IHC

**Recommended Dilution:** Flow Cytometry: < / = 1  $\mu$ g/10e6 cells.

Immunohistochemistry (Acetone-Fixed, Frozen Tissue Sections).

Reactivity:MouseHost:MouseIsotype:IgG2a

Clonality: Monoclonal

Immunogen: C3HalphaBDF1 mouse splenocytes.

**Specificity:** This antibody is specific to an epitope in the alpha 3 domain that is common to H-2Kd and H-

2Dd. It binds to a common determinant in the alpha 3 domains of H-2Kd and H-2Dd in the

presence or absence of beta 2 Microglobulin. (Ref.7,8) It cross reacts with the alpha 3 domain of H-2Kb. (Ref.1,2)

**Formulation:** PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.

Label: PE

State: Liquid purified Ig fraction.

Label: R-Phycoerythrin

**Concentration:** lot specific

**Purification:** Liquid purified lg fraction.

Conjugation: PE

**Storage:** Store the antibody undiluted at 2-8°C.

DO NOT FREEZE!

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.







## Background:

The 'classical' MHC Class I molecules are histocompatibility antigens encoded by the H-2 gene complex and consist of heterodimers of highly polymorphic alpha chains noncovalently associated with the invariant beta 2-Microglobulin. (Ref.3,4) These antigens are expressed on most nucleated cells but expression varies on different cell types. MHC Class I molecules present endogenously synthesized peptides to CD8+ T lymphocytes, which are usually cytotoxic T cells. (Ref.5) MHC Class I antigens expressed on thymic epithelial cells regulate the positive and negative selection of CD8+ T cells during T cell ontogeny. (Ref.3,6)