

## Product datasheet for AM08084LE-N

## Cd86 Rat Monoclonal Antibody [Clone ID: 2D10]

## **Product data:**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Clone Name:	2D10
Applications:	FC, FN, IHC, IP
Recommended Dilution:	<ul> <li>ELISA. (Ref.4).</li> <li>Flow Cytometry: Identification and enumeration of CD86+ cells. (Ref.1,2,4,6,7)</li> <li>Immunoprecipitation (Ref.1)</li> <li>Immunohistochemistry on Frozen Sections: Identification of CD86+ cells. (Ref.4).</li> <li>In vivo and in vitro blocking assays. (Ref.1-3,5)</li> </ul>
Reactivity:	Mouse
Host:	Rat
lsotype:	lgG2b
Clonality:	Monoclonal
Immunogen:	Mouse B cell lymphoma cell line 5C2
Specificity:	This antibody recognises CD86, the B7-2 costimulatory molecule and a ligand for CD28 and CD152.
Formulation:	PBS without preservatives. State: Low Endotoxin State: Liquid purified lg fraction. Stabilizer: None Preservative: None
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD86 antigen
Database Link:	Entrez Gene 12524 Mouse P42082



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Cd86 Rat Monoclonal Antibody [Clone ID: 2D10] – AM08084LE-N
Background:	CD86, also known as B72, is a type I transmembrane glycoprotein and a member of the immunoglobulin superfamily of cell surface receptors. It is expressed at high levels on resting peripheral monocytes and dendritic cells and at very low density on resting B and T lymphocytes. CD86 expression is rapidly upregulated by B cell specific stimuli with peak expression at 18 to 42 hours after stimulation. CD86, along with CD80/B71, is an important accessory molecule in T cell costimulation via its interaction with CD28 and CD152/CTLA4. Since CD86 has rapid kinetics of induction, it is believed to be the major CD28 ligand expressed early in the immune response. It is also found on malignant Hodgkin and Reed Sternberg (HRS) cells in Hodgkin's disease. CD86 interacts with HHV8 (KSHV) MIR2 protein.
Synonyms:	CD28LG2, Activation B7-2 antigen, B70, BU63, B7.2, FUN-1

~ \$%~

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US