

Product datasheet for **AM08088FC-N**

Lamp2 Rat Monoclonal Antibody [Clone ID: ABL-93]

Product data:

Product Type:	Primary Antibodies
Clone Name:	ABL-93
Applications:	FC
Recommended Dilution:	Flow Cytometry: < / = 1 µg/10e6 cells.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Glycoproteins purified from BALB/c mouse embryo 3T3 cell line. (Ref.1)
Specificity:	This antibody is specific to CD107b/LAMP-2 (Mr 110-120 kDa with a core protein of Mr ~40 kDa).
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	lysosomal-associated membrane protein 2
Database Link:	Entrez Gene 16784 Mouse P17047



[View online »](#)

Background:

CD107b, also known as lysosomal-associated membrane protein 2 (LAMP-2), is a heavily glycosylated, type I transmembrane protein that constitutes one of the two major sialoglycoproteins on lysosomal membranes that can be used to distinguish lysosomes from other organelles. (Ref.1) CD107b is similar in structure to Mac-1 and is expressed primarily on epithelial cells and macrophages in a variety of organs in normal and Beige mutant mice. (Ref.2)

Synonyms:

LAMP-2, LAMP-2C, LAMPB