

## Product datasheet for AM08090BT-N

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 7E7]

**Product data:** 

**Product Type:** Secondary Antibodies

Product Name: Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 7E7]

Clone Name: G1 7E7
Applications: FC, IF

**Recommended Dilution:** Flow Cytometry:  $\leq 1 \mu g/106$  cells.

Immunofluorescence Microscopy: Identification and enumeration of IgG1+ cells.

**ELISA:** 1/5,000-1/20,000.

Reactivity: Rat

Host: Mouse

Immunogen: Rat IgG1 hybridoma

Isotype: IgG1

**Formulation:** PBS containing 0.09% Sodium Azide as preservative.

Label: Biotin

State: Liquid purified Ig fraction.

**Concentration:** lot specific

Conjugation: Biotin

Storage: Store undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.