

Product datasheet for AM08090HR-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 7E7]

Product data:

Product Type: Secondary Antibodies

Product Name: Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 7E7]

Clone Name: G1 7E Applications: ELISA

Recommended Dilution: ELISA: 1/2,000-1/5,000.

Reactivity: Rat
Host: Mouse

Isotype: IgG1

Formulation: 50% Glycerol/50% PBS, pH 7.4 without preservatives.

Label: HRP

State: Liquid purified Ig fraction. Label: Horseradish Peroxidase

Conjugation: HRP

Storage: Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

