

## Product datasheet for AM08091AP-N

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 1C5]

**Product data:** 

**Product Type:** Secondary Antibodies

Product Name: Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 1C5]

Clone Name: G1 1C Applications: ELISA

Recommended Dilution: ELISA: 1/2,000-1/4,000.

**Reactivity:** Rat

Host: Mouse Isotype: IgG2b

**Formulation:** 50mM Tris/1mM MgCl2/50% Glycerol, pH 8.0 containing 0.09% Sodium Azide as preservative.

Label: AP

State: Liquid purified lg fraction. Label: Alkaline Phosphatase

Conjugation: AF

Storage: Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

