

## Product datasheet for **AM08091C5-N**

### Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 1C5]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 1C5]
Clone Name:	G1 1C5
Applications:	FC, IF
Recommended Dilution:	<b>Flow Cytometry:</b> < / = 0.3 µg/10e6 cells. <b>Immunofluorescence Microscopy:</b> Identification and enumeration of IgG1+ cells.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2b
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Cy5 State: Liquid purified Ig fraction. Label: Cyanine 5
Concentration:	lot specific
Conjugation:	Cy5
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)