

Product datasheet for **AM08091FC-N**

Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 1C5]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 1C5]
Clone Name:	G1 1C5
Applications:	FC, IF
Recommended Dilution:	FLISA: 1/200-1/400. Flow Cytometry: 1 ≤ µg/10 ⁶ cells. Immunofluorescence Microscopy: Identification and enumeration of IgG1+ cells.
Reactivity:	Rat
Host:	Mouse
Immunogen:	Rat IgG1 hybridoma
Isotype:	IgG2b
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)