

Product datasheet for **AM08091RP-N**

Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 1C5]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: G1 1C5]
Clone Name:	G1 1C5
Applications:	FC, IF
Recommended Dilution:	Flow Cytometry: < / = 0.02 µg/10e6 cells. Immunofluorescence Microscopy: Identification and enumeration of IgG1+ cells.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2b
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)