

## Product datasheet for **AM08092AP-N**

### Rat IgG2a (gamma 2a chain specific) Mouse Monoclonal Antibody [Clone ID: 2A 8F4]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2a (gamma 2a chain specific) Mouse Monoclonal Antibody [Clone ID: 2A 8F4]
Clone Name:	2A 8F4
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/2,000-1/4,000.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Formulation:	50mM Tris/1mM MgCl <sub>2</sub> /50% Glycerol, pH 8.0 containing 0.09% Sodium Azide as preservative. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)