

## Product datasheet for **AM08092BT-N**

### Rat IgG2a (gamma 2a chain specific) Mouse Monoclonal Antibody [Clone ID: 2A 8F4]

#### Product data:

|                       |   |
|-----------------------|---|
| Product Type:         | Secondary Antibodies  |
| Product Name:         | Rat IgG2a (gamma 2a chain specific) Mouse Monoclonal Antibody [Clone ID: 2A 8F4]  |
| Clone Name:           | 2A 8F4  |
| Applications:         | ELISA, FC, IF   |
| Recommended Dilution: | <b>Flow Cytometry:</b> $\leq 1 \mu\text{g}/10^6$ cells.<br><b>Immunofluorescence Microscopy:</b> Identification and enumeration of IgG2a+ cells.<br><b>ELISA:</b> 1/5,000-1/20,000. |
| Reactivity:           | Rat   |
| Host:                 | Mouse   |
| Immunogen:            | Rat IgG2a hybridoma.  |
| Isotype:              | IgG1  |
| Formulation:          | PBS containing 0.09% Sodium Azide as preservative.<br>Label: Biotin<br>State: Liquid purified Ig fraction.  |
| Concentration:        | lot specific  |
| Conjugation:          | Biotin  |
| Storage:              | Store undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing.  |



[View online »](#)