

Product datasheet for **AM08092FC-N**

Rat IgG2a (gamma 2a chain specific) Mouse Monoclonal Antibody [Clone ID: 2A 8F4]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2a (gamma 2a chain specific) Mouse Monoclonal Antibody [Clone ID: 2A 8F4]
Clone Name:	2A 8F4
Applications:	FC, IF
Recommend Dilution:	Flow Cytometry: < / = 1 µg/10e6 cells. Immunofluorescence Microscopy: Identification and enumeration of IgG2a+ cells.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)