

## Product datasheet for **AM08093HR-N**

### Rat IgG2b (Gamma-2b chain specific) Mouse Monoclonal Antibody [Clone ID: 2B 10A8]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2b (Gamma-2b chain specific) Mouse Monoclonal Antibody [Clone ID: 2B 10A8]
Clone Name:	2B 10A8
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/4,000-1/8,000.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2b
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)