

## Product datasheet for **AM08093PU-N**

### Rat IgG2b (Gamma-2b chain specific) Mouse Monoclonal Antibody [Clone ID: 2B 10A8]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2b (Gamma-2b chain specific) Mouse Monoclonal Antibody [Clone ID: 2B 10A8]
Clone Name:	2B 10A8
Applications:	ELISA, FC, IF
Recommended Dilution:	<b>Flow Cytometry.</b> <b>Immunofluorescence Microscopy:</b> Identification and enumeration of IgG2b+ cells. <b>ELISA.</b> <b>Second step reagent</b> for rat IgG2b Monoclonal Antibodies.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2b
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)