

## Product datasheet for **AM08093RP-N**

### Rat IgG2b (Gamma-2b chain specific) Mouse Monoclonal Antibody [Clone ID: 2B 10A8]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2b (Gamma-2b chain specific) Mouse Monoclonal Antibody [Clone ID: 2B 10A8]
Clone Name:	2B 10A8
Applications:	FC, IF
Recommended Dilution:	<b>Flow Cytometry:</b> < / = 0.02 µg/10e6 cells. <b>Immunofluorescence Microscopy:</b> Identification and enumeration of IgG2b+ cells.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2b
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)