

Product datasheet for **AM08094BT-N**

Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]
Clone Name:	2C 8F1
Applications:	ELISA, FC, IF
Recommended Dilution:	Flow Cytometry: $\leq 1 \mu\text{g}/10^6$ cells. Immunofluorescence Microscopy: Identification and enumeration of IgG2c ⁺ cells. ELISA: 1/5,000-1/10,000.
Reactivity:	Rat
Host:	Mouse
Immunogen:	Rat IgG2c hybridoma.
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative Label: Biotin State: Liquid purified Ig fraction
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)