

Product datasheet for **AM08094C5-N**

Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]
Clone Name:	2C 8F1
Applications:	FC, IF
Recommended Dilution:	Flow Cytometry: < / = 0.3 µg/10e6 cells. Immunofluorescence Microscopy: Identification and enumeration of IgG2c+ cells.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Cy5 State: Liquid purified Ig fraction. Label: Cyanine 5
Concentration:	lot specific
Conjugation:	Cy5
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)