

Product datasheet for **AM08094HR-N**

Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]
Clone Name:	2C 8F1
Applications:	ELISA
Recommended Dilution:	ELISA: 1/4,000-1/8,000.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2a
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)