

## Product datasheet for **AM08094PU-N**

### Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]
Clone Name:	2C 8F1
Applications:	ELISA, FC, IF
Recommended Dilution:	<b>Flow Cytometry.</b> <b>Immunofluorescence Microscopy:</b> Identification and enumeration of IgG2c+ cells. <b>ELISA.</b> <b>Second step reagent</b> for rat IgG2c Monoclonal Antibodies.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2a
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)