

Product datasheet for **AM08094RP-N**

Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG2c (gamma 2c chain specific) Mouse Monoclonal Antibody [Clone ID: 2C 8F1]
Clone Name:	2C 8F1
Applications:	ELISA, FC, IF
Recommended Dilution:	Flow Cytometry: FITC and BIOT conjugates $\leq 1 \mu\text{g}/10^6$ cells. Immunofluorescence Microscopy: Identification and enumeration of IgG2c+ cells. ELISA.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)