

Product datasheet for **AM08095BT-N**

Rat IgM (Mu chain specific) Mouse Monoclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgM (Mu chain specific) Mouse Monoclonal Antibody
Applications:	ELISA, FC, IF
Recommended Dilution:	Flow Cytometry (< / = 1.0 µg/10e6 cells). Immunofluorescence Microscopy: Identification and enumeration of IgM+ cells. ELISA. Second step reagent for rat IgM Monoclonal Antibodies. Surface IgM Detection. B Cell Stimulation.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C.



[View online »](#)