

Product datasheet for AM08095C5-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: M 2A1]

Product data:

Product Type: Secondary Antibodies

Product Name: Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: M 2A1]

Clone Name: M 2A1
Applications: FC, IF

Recommended Dilution: Flow Cytometry (< $/ = 0.3 \mu g/10e6$ cells) and Immunofluorescence Microscopy:

Identification and enumeration of IgM+ cells.

Reactivity: Rat

Host: Mouse Isotype: IgG1

Formulation: PBS containing 0.09% Sodium Azide as preservative.

Label: Cy5

State: Liquid purified Ig fraction.

Label: Cyanine 5

Concentration: lot specific

Conjugation: Cy5

Storage: Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.