

## Product datasheet for AM08095FC-N

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: M 2A1]

**Product data:** 

**Product Type:** Secondary Antibodies

Product Name: Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: M 2A1]

Clone Name: M 2A1
Applications: FC, IF

**Recommended Dilution:** Flow Cytometry: < / = 1.0  $\mu$ g/10e6 cells.

**Immunofluorescence Microscopy:** Identification and enumeration of IgM+ cells.

**Reactivity:** Rat

Host: Mouse

**Immunogen:** Rat IgM hybridoma.

Isotype: IgG1

**Formulation:** PBS containing 0.09% Sodium Azide as preservative.

Label: FITC

State: Liquid purified Ig fraction.

Label: Fluorescein Isothiocyanate Isomer 1

**Concentration:** lot specific

Conjugation: FITC

Storage: Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.

