

Product datasheet for **AM08095HR-N**

Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: M 2A1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: M 2A1]
Clone Name:	M 2A1
Applications:	ELISA
Recommended Dilution:	ELISA: 1/2,000-1/5,000.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)