

Product datasheet for **AM08095PU-N**

Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: M 2A1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: M 2A1]
Clone Name:	M 2A1
Applications:	ELISA, FC, IF
Recommended Dilution:	Flow Cytometry and Immunofluorescence Microscopy: Identification and enumeration of IgM+ cells. ELISA. Second step reagent for rat IgM Monoclonal Antibodies. Surface IgM Detection. B Cell Stimulation.
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG1
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)